

Key Features

- Precision ECG R-wave peak detection
- Synchronized ECG trigger outputs
- 4 lead ECG configuration with auto lead selection
- Built-in ECG simulator
- On-screen color coded trigger pulse indication
- Ethernet communication interface
- X-Ray on interface
- Analog ECG waveform output
- Digital ECG waveform output & storage
- Multi-Language user interface
- Optional strip chart recorder
- Patient isolation/protection
- Universal power supply/voltage
- FDA 510(k) & CE Mark

Product Description


The Ivy Biomedical Systems' Model 7800 is our fifth generation of cardiac gating monitors. It is ideal for use in applications requiring precision ECG R-wave synchronization, such as computed tomography (CT) imaging systems for cardiac and vascular studies.

Value added features include: auto ECG lead selection ensures the best trigger vector will be used; built-in ECG simulator for pre-scan testing of the entire system; electrode impedance check for reliable scans; X-Ray on signal input enables ECG waveform file capture for post scan analysis; two channel, electrically isolated Ethernet port allows for data interface to both the console and gantry display. An optional strip chart recorder is also available for hardcopy documentation of ECG rhythms.

Ivy Biomedical System's renowned high quality and reliability ensure consistent uptime and long operational life.

Synchronized ECG Trigger Output*

Trigger Delay*	< 2ms*
R-to-R Accuracy	±75µs dither (typ.) @ 1mV input
Pulse Width**	1, 50, 100 or 150 ms
Pulse Amplitude**	0V to +5V or -10V to +10V
Pulse Polarity**	Positive or Negative
Source Current:	8mA @ 3V (Max)
Sink Current:	16mA (Max)



* Input signal test conditions: ½ sine wave, 60ms width, 1mV amplitude, 1 pulse/sec; Clinical Settings: 2-8ms delay dependant on QRS width

** Pre-configured at the factory only

ECG

Configuration	4-Lead system
Trigger Lead Selection	I, II, III, or AUTO
Second Lead Display	I, II, III
ECG Simulator	Integrated
Patient Isolation	>4 kV rms, 5.5 kV peak
Frequency Response	0.67 - 100Hz unfiltered 1.5-40Hz Filtered
Notch Filter	50/60 Hz (auto)
CMRR	≥ 90dB
Tall T-wave Rejection	≤ 1.2 * R-wave
Pacer Rejection (user on/off)	0.1 to 2ms pulse width @ ±2 to ±700 mV
Defibrillator Protection	360 J discharge; < 5 sec recovery time (Type CF)

Cardiotach

Adult	10-300 bpm
Pediatric/Neonate	10-350 bpm
Accuracy	±1% ±1 bpm
Resolution	1 bpm
Sensitivity	300 µV peak
HR Averaging	Exponential @ 1Hz; 2 or 8 sec max response time

Alarms

High HR Limit	15-250 bpm (5 bpm inc.)
Low HR Limit	10-245 bpm (5 bpm inc.)
Asystole	R-to-R interval > 6 sec
ECG Lead Off	Each detached lead
Check ECG Lead	Lead imbalance > 0.5V

Manufactured by:



Ivy Biomedical Systems, Inc.
11 Business Park Drive
Branford, Connecticut 06405 USA
Toll Free 800 247 4614
Main 203 481 4183
Fax 203 481 8734
www.ivybiomedical.com

Display

Waveform Type	Dual trace; Freeze Active Matrix TFT Color Touch Screen LCD
Resolution Size	640x480 pixels 8.4" (21.3 cm) diagonal

Input/Output Interface

Synch Output	BNC; Provides trigger pulse output synch to ECG R- wave peak
ECG Output	¼" stereo jack; Provides trigger pulse output synchro- nized to ECG R-wave peak as well as analog ECG waveform output
Aux Conn	DB-9; X-Ray On interface; ECG trigger pulse; ECG waveform
RS-232 Comm	Micro DB-9; device interface
Ethernet Comm	RJ-45; 10BaseT, IEEE 802.3 (2 channels)
Data Storage	USB 2.0 Type A 200 ECG Events (FIFO)

Mechanical

Size (HxWxD)	22.1x23.5x15.5 cm (8.7x9.25x6.1 inches)
Weight	2.54 kg (5.6 lbs.)
Case Material	Polycarbonate

Electrical

Input Voltage	100-120Vac; 200-230Vac
Frequency	50/60 Hz
Power Consumption	45 VA (max.)
Power Recovery	Auto if power restored within 30 seconds

Environmental

Water Resistance	IPX1
Operating	
Temperature Range	5°C to 40°C
Relative Humidity	0% to 90% non-condensing
Altitude	-100m to +3,600m
Storage	
Temperature Range	-40°C to +70°C
Relative Humidity	5% to 95% non-condensing
Altitude	-100m to +14,000m

Options

Integrated Recorder	2 trace, direct thermal
Mounting Plate	3" adaptor for rollstand with 3" receiver plate
Roll Stand	with 3" receiver plate

(Specifications subject to change without notice)

Distributed by:

Accessories

Electrodes	Low impedance; 10% KCl wet gel sponge type
ECG Leads	4-lead carbon (radio translucent) with pinch clips; AHA or IEC color code; 30" or 36" lengths available
Trunk Cable	40", 5', 10' cable with 6-pin AAMI connector

Globalization

User Interface	12 selectable languages
Operator's Manual	33 languages on CD
Registrations	Multiple countries

Compliance & Certifications

ANSI/AAMI ES60601-1: A1:2012, C1:2009/
(R)2012 and A2:2010/(R)2012
IEC 60601-1 Edition 3.1 (2012)/EN 60601-1:2006
+ A1:2013 + A12:2014
IEC 60601-1-2:2014+AMD1:2020, Edition 4.1
IEC 60601-1-6:2010 (Third Edition) + A1:2013; IEC
62366:2007 (First Edition) + A1:2014
IEC 60601-1-8:2006 (Second Edition) + Am.1:2012
IEC 60601-2-27 (2011)
IEC 62304:2006
CAN/CSA-C22.2 No. 60601-1:2014
CAN/CSA-C22.2 No. 60601-1-2:2016
EU MDR 2017/745
CE 2862
ISO 13485:2016
RoHS 2011/65/EU
WEEE 2012/19/EU
FDA/CGMP
MDSAP

Notified Body

Intertek Medical Notified Body AB,
Identification Number 2862
MDR Classification IIb

Authorized Representative

EC REP: Emergo Europe
CH REP: MDSS CH GmbH



For additional specifications,
refer to Operator Manual

